

Instructions for assembly high brakes extended moved ahead

How to assemble the extended moved ahead high brakes:

1. When the brakes are delivered to you they are assembled as shown in figure 1. Before you start to assemble the brakes you should check that the spring (3) is not trapped between the clamp (4) and the spacer (7). If it is trapped, you should reposition the spring round the spacer.




Figure 1.

2. Loosen the bolt slightly with the socket head (6) (figure 2) to enable you to place the brake clamps (4 and 5) round the seat tubing on the wheelchair. (See figure 3).

3. Then try tightening the bolt (6) and move the brake handle (1) forward to brake. The brakes should be adjusted so that they sink into the tyres by around 5mm when in locked position.

- If you feel that the brakes are too loose, push the brake handle backwards and loosen the bolt (6) moving the brake back slightly. Screw in the bolt securely and check that you are satisfied with the brakes.

- If you feel the brakes are too stiff, loosen the bolt (6) and move the brake forwards slightly. Screw in the bolt securely and check that you are satisfied with the brakes.

4. Make sure that the brake locking block (2) is positioned as straight as possible against the rear wheel for the best possible braking effect.




Figure 2.

5. Repeat the same procedure on the other side.




Figure 3.