

2.1 KILOGRAMS





# PANTHERA X

PANTHERA X™ is by far the lightest and most high-tech wheelchair ever made. This chair is built from the most sophisticated materials using the most advanced technology available today.

The result is a chair with extremely low weight, but still without compromising stability and comfort. As users ourselves, we can definitely say that PANTHERA X™ is a dream come true for active users.



The PANTHERA X™  
you see here weighs  
just 4.4 kilos.

# WEIGHT

## 2.1 KILOS

Lifting PANTHERA X™ into your car is unreal. With a transport weight of 2.1 kilos, it feels like lifting paper.

## 4.4 KILOS

With a total weight of 4.4 kilos, PANTHERA X™ is amazingly light. There's no other chair that even comes close.



PANTHERA X TECHNOLOGY

PANTHERA X™  
is based on the new  
X technology.

## STABILITY

A rigid frame is the key to easy handling. All flex requires energy, energy that comes from the user. By using stiff materials and a box frame design, PANTHERA X™ reduces energy losses to nearly zero. PANTHERA X™ handles and rolls so effortlessly, you just have to drive it to believe it.

## COMFORT

Good seat position, especially proper back support, are vital to the everyday user. On PANTHERA X™, you can adjust the backrest angle as well as fine-tune the backrest upholstery.

PANTHERA  TECHNOLOGY

• The real challenge behind  
• PANTHERA X™ has been  
• determining the right combination  
• of fiber qualities and fiber directions.

• PANTHERA X™ is just  
• unbelievably strong.



# MATERIALS

As you can see PANTHERA X™ is made out of carbon fiber but beneath the surface lies its secret. Making something out of carbon fiber is difficult, but making it extremely light, yet at the same time strong, is the real challenge.

It has taken us five years to develop PANTHERA X™. We've spent endless hours of testing, just to optimize the fiber directions and determine the right combination of fiber qualities.

## FRAME

Carbon fiber

## BACKREST

Carbon fiber

## FOOTREST

Aluminum 7075

## PUSHRIM

Aluminum 7075




# PARTS

On PANTHERA X™ we've basically made every part by ourselves. Even the smallest parts, like bolts and nuts, have been redesigned to save even the smallest fraction of a gram.

Time, money or effort has never been an issue.

That's PANTHERA X™ TECHNOLOGY



• Almost everything you see on  
• PANTHERA X™ we've  
• made ourselves.

[WWW.PANTHERA-X.COM](http://WWW.PANTHERA-X.COM)

**panthera.**

[WWW.PANTHERA.SE](http://WWW.PANTHERA.SE)

Panthera AB, Gunnebogatan 26, SE-163 53 SPÅNGA, Tel: +46 8 761 50 40, Fax: +46 8 621 08 10, [panthera@panthera.se](mailto:panthera@panthera.se)