

Adjusted by:

E-mail:

Tel:

Date:

To mount the Mototronik on a Panthera wheelchair use the Quick Connector attachment. Select your wheelchair model and seat width to get the right part number for your Quick Connector when ordering.

Select your Quick Connector \*- See table to the right

- ST Quick Connector Mototronik X L1
- ST Quick Connector Mototronik X L2
- ST Quick Connector Mototronik X L3
- ST Quick Connector Mototronik 19 L1
- ST Quick Connector Mototronik 19 L2
- ST Quick Connector Mototronik 19 L3
- ST Quick Connector Mototronik 19 L4

CLEAR

\* Works only with Mototronik 12"

	Panthera X	U2 / U3 Light	S3
<b>X L1</b>	33 cm, 36 cm		
<b>X L2</b>	39 cm, 42 cm		
<b>X L3</b>	45 cm		
<b>19 L1</b>		33 cm, 36 cm	
<b>19 L2</b>		39 cm, 42 cm	33 cm
<b>19 L3</b>		45 cm	36 cm, 39 cm
<b>19 L4</b>			42 cm, 45 cm

Art.nr:

