

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

S3 0° has a ajustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 has a seat height of 45 cm in the front and 45 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options. **Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45-50 cm.**

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. **Blue, underlined text is clickable to show info** or image. **BB** = seatwidth, **HH**= backrest height

**Original design**

**Alternative options**

**Other options**

[info](#)



**1) Chassis-colour**

Click on the link to see up-to-date colours.

Fill in art.no. in the square

**2) Seat width (BB)**

[art.no.](#)

Clear

- 33 cm  36 cm  39 cm  42 cm  45 cm  50 cm

Art.no include: chassis, seat, calf strap, backrest lock

**3) Seat depth**

[art.no.](#)

Clear

- 35 cm  40 cm

- 37,5 cm  42,5 cm

**4) Backrest**

[Standard](#)  
523BBHH VER2

Clear

Height (HH)

- 20 cm  35 cm  
 25 cm  40 cm  
 30 cm  45 cm  
 Clear

[Tapered](#)

525BBHH VER2

[Curved back](#)

524BBHH VER2

[Backrest, extended range](#)

[Frame for fixed backrest systems](#)

art. \_\_\_\_\_

**5) Rear axle**

[Carbon standard](#)

33023 BB

Clear

[Carbon +15 mm](#)

33023BC\*

\* seat width **BB** +1

Use this rear axle when chosing arm support at Accessories and/or increased seat height

[Other](#) rear axle or seat height

[Increase the seat height](#)

**6) Castor and fork**

[120 mm](#)

Clear

[90 mm](#)

[150 mm](#)

**7) Footrest**

[Footrest Titanium](#)

36001BB

Clear

[Extended, reinforced](#)

36111BB

[Foldable footplates](#)

[art.no.](#)

[Heel strap for footplates](#)

[art.no.](#)

[Plastic plate for footrest](#)

36010BB

[Other footrest](#)

**8) Brake**

[S3 standard](#)

4630000

Clear

[One arm brake, right](#)

[art.no.](#)

[One arm brake, left](#)

[art.no.](#)

**9) Rear wheel - without pushrim**

[24" standard](#)

[art.no.](#)

Clear

*Rear wheel sizes*

24" [Spinergy X](#) black spokes

24" [Spinergy X](#) white spokes

25" [Spinergy X](#) black spokes

25" [Spinergy X](#) white spokes

26" [Spinergy X](#) black spokes

26" [Spinergy X](#) white spokes

[art.no.](#)

[Other rear wheel](#)

**Tire - choose tire**

[Schwalbe Right Run](#) 24",25",26"

Clear

[Rugged](#) 24"

Schwalbe [Marathon Plus](#) 24",25",26"

Schwalbe [One](#) 24",25"

[PU](#) 24" (massive)

[reTyre](#)-tyre skins (does not fit PU tyres)

**10) Pushrim**

[Titanium](#)

190DD11 (x2)

Clear

**DD** = Wheel diameter

[Alu X](#) 24"

1904404 (x2)

[Tetragrip](#)

19052DDBT x2

[Paragrip](#)

19051DDBT x2

[Alu X](#) 25"

1904405 (x2)

[Alu X](#) 26"

1904406 (x2)

**DD** = Wheel diameter

[Maxgrip](#) 24" (contains latex)

1903060 x2

[Maxgrip](#) 22" (contains latex)

1903080 x2

## Accessories



### Side guard

- [Side guard S3](#) [Video](#)  
[art.no.](#)
- Clear



### Side guard Carbon

- [Carbon](#)  
[art.no.](#)
  - [Carbon, small cover](#)  
[art.no.](#)
  - [Carbon, large cover \\*](#)  
[art.no.](#)
- \* Only with one arm brake



### Arm support

- [Arm support S3](#) [Video](#)  
[art.no.](#)



### Push handles H/S

- [Height adjustable S3](#) [Video](#)  
4252060 VER2
- [Height adjustable S3+10 cm](#)  
4252063 VER2
- Max. userweight 100kg
- Clear



### Push handle Foldable

- [Foldable S3](#) [Video](#)  
4245000



### Push bar

- [Push bar S3](#)  
42400BB VER2
- Max. userweight 100kg



### Anti tip

- [Anti tip](#) std seat height  
4353080 [Video](#)

[Other seat height](#)



### Cushion

Seat depth 35cm    Seat depth 40cm

- [2,5 cm](#)     [2,5 cm](#)
- 572BB35     572BB40
- [5 cm](#)     [5 cm](#)
- 575BB35     575BB40
- Clear



### Spoke protectors

- [24" logo](#)  
4500000
- [24" no logo](#)  
4500011
- Clear



### Other

- [Tube protection](#) \_\_\_\_\_ p  
9900251 (x3/wheelchair)
- [Crutch holder S3](#) \_\_\_\_\_ p  
9900301 (x1/crutch) [Video](#)
- [Backrest wedge](#) \_\_\_\_\_ p  
4710012

Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)

- [Freewheel](#)



### Wedge cushion

- [Wedge cushion](#)  
57710BB

## Set-up adjustments



### Back rest angle

Note measure (A) for desired back rest angle.  
14-28 mm. Standard=20 mm

A= \_\_\_\_\_ mm



### Rear axle adjustment

Note measure (A) for desired balance.  
Standard=60 mm


A= \_\_\_\_\_ mm



### Footrest height [Footbow- Footplates](#)

Note measure (A) and (B) for desired footrest.  
[Bambino foot plate with Feewheel - info](#)

A= \_\_\_\_\_ mm    B= \_\_\_\_\_ mm



### Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle

Art.no.9000152

NOTE! See [manual](#) for more information.

Extra