

Panthera U3

Framgaffel med länkhjul vid olika storlekar på drivhjul

90 mm

Drivhjul	Framgaffel	art.nr.		
24"		2100030	Gaffel 90, hjul 90, tapp 3	x2
25"		2100060	Gaffel 105, hjul 90, tapp 3	x2
26"		2100090	Gaffel 120, hjul 90, tapp 3	x2

120 mm

Drivhjul	Framgaffel	art.nr.		
24"		2100020	Gaffel 75, hjul 120, tapp 3	x2
25"		2100040	Gaffel 90, hjul 120, tapp 3	x2
26"		2100070	Gaffel 105, hjul 120, tapp 3	x2

ENG

Castor fork with castor for different rearwheel sizes.

90 mm

Rearwheel	Castor fork	art.no.		
24"		2100030	Fork 90 Castor 90 Bolt 3	x2
25"		2100060	Fork 105 Castor 90 Bolt 3	x2
26"		2100090	Fork 120 Castor 90 Bolt 3	x2

120 mm

Rearwheel	Castor fork	art.no.		
24"		2100020	Fork 75 Castor 120 Bolt 3	x2
25"		2100040	Fork 90 Castor 120 Bolt 3	x2
26"		2100070	Fork 105 Castor 120 Bolt 3	x2

120 mm