penthere.

S3 series



A NEW STANDARD

At Panthera, we spend all day, every day, making life easier for wheelchair users.

Now, we proudly present the S3 series. A standard chair that is everything but.

113 IMPROVEMENTS TO EVERYDAY LIFE

Our mindset, as always, has been to improve, simplify and make possible.

Carbon fibre, chromemolybden and titanium provide rigidity, strength and low weight. Integrated functionality and fewer parts make way for easy assembly and adjustment. The driving is smoother, the lifting lighter and the sitting more comfortable. In fact, we found 113 ways to improve everyday life. Enjoy the ride.

Accessories and functions



SIDE GUARDS Soft upper part creates new possibilities



ARM SUPPPORT
Fully adjustable by the user without tools



ANTI TIP DEVICE
Activated and deactivated by the user sitting in the wheelchair





CASTOR WHEELS
Unique design combines two functions in one.
Easy rolling and also stability on soft surface.



ADJUSTABLE PUSH HANDLES Fully adjustable without tools



FOLDABLE PUSH HANDLES



penthere.

Models



Panthera S3 S3, S3 short, S3 short-low, S3 large



Panthera U3





Panthera S3 Swing S3 Swing, S3 Swing-short



Panthera U3 Light / L



facts	S3	U3	SWING	U3 LIGHT
seat width	33-50 cm	33-45 cm	33-50 cm	33-45 cm
	11-20"	13-18"	13-20"	13-18"
total weight	8.3 kg*	8.2 kg*	9.4 kg*	6.1 kg*
	18 lb	18 lb	18 lb	13 lb
transport weight	4.7 kg*	4.6 kg*	4.6 kg**	3.8 kg*
	10 lb	10 lb	10 lb	8 lb
max. user weight	150 kg	150 kg	150 kg	100 kg
	330 lb	330 lb	330 lb	220 lb

*with brake

**without leg rests

All models are tested and approved according to standard EN12183:2014 and ISO-7176-16 S3, U3 and Swing are tested and approved according to crash test standard ISO-7176-19:2001

