

NEW ATTACHMENT FOR PUSHRIM MAXGRIP 22"

Click on images for larger format



OLDER ATTACHMENT



NEW ATTACHMENT

New attachment for pushrim MaxGrip 22"

New design of the attachment plates on Pushrim MaxGrip 22".

We introduce this improvement to get the same distance between the rim and the pushrim on Pushrim MaxGrip 22" as all other pushrims in our range.

Art.no.	Description
1903080	Pushrim MaxGrip 22" wheel 6P nuts

Contact Panthera if you have questions.

Contact person at Panthera:

Export / Marketing Manager

Helena Lemoine

Tel: +46-557 608 77

E-mail: helena.lemoine@panthera.se