

Phasing out "High brake moved ahead" for Bambino 1

Click on images for larger format



HIGH BRAKE MOVED AHEAD

Phasing out "High brake moved ahead" for Bambino 1

We will now start phasing out "High brake moved ahead", used for Bambino 1 and replace it with "Brake Bambino 3".

Valid for brake: [High brake moved ahead](#)

It will be replaced with: [Brake Bambino 3](#), [Brake Bambino 3 moved forward](#) and [Adjustment plate Bambino 3 moved backwards](#) depending on rearwheel size and the position of the rear axle, center of gravity (COG).

The price for "Brake Bambino 3" is lower than for "High brake moved ahead".



BRAKE BAMBINO 3



BRAKE BAMBINO 3
MOVED FORWARD



ADJUSTMENT PLATE BAM-
BINO 3 MOVED BACKWARDS

Phasing out		Replaced with	
4600050	High brakes pair moved ahead	4630200	Brake B3 (B3 =Bambino 3) cpl pair
4600051	High brakes left moved ahead	4630201	Brake B3 cpl left
4600052	High brakes right moved ahead	4630202	Brake B3 cpl right
		4630300	Brake B3 moved forward cpl pair
		4630301	Brake B3 moved forward cpl left
		4630302	Brake B3 moved forward cpl right
		4630236	Adjustment plate Brake B3 pair
		4630235	Adjustment plate Brake B3