

# Info for Tray table Bambino 3

Click on images for larger format



TRAY TABLE BAMBINO 3

## Info for Tray table Bambino 3

Panthera are working on an improvement regarding the attachment of the Tray Table.

From today, if ordering Tray Table cpl as after market accessory or spare part Kit 2 you need to state what seat width wheelchair the Tray will be used on.

Once the improvement is implemented we will inform you of this.

Brickbord / Tray

Magnetfäste / Magnet holder

Rev. 2020-01

POS	DESCRIPTION	ART.NO	QTY
	Tray table kit 1	6810015	KIT
	Tray table kit 2	6810016	KIT
	Holder magnetic traytable B3 cpl	6810010	KIT
1	Tray table B3	6810001	1
2	Stålkapp Brickbord fzb Steelcup, tray fzb	6810022	3
3	Plugg Brickbord Plug, tray	6810003	3
4	Brace holder tray	6810004	2
5	Busching Tray	6810005	2
6	Brace traytable	6810006	2
7	O-ring Tray	6810012	3
8	Lock mit DIN 985 steel M5 fzb	7000350	4
9	Bolt MLCBS M5x10 fzb	6002508	4
10	Magnet holder	6810008	1
11	Magnet	9900307	3

Click on images for larger format

Please do not hesitate with any questions:

Commercial Support Sales

**Therese Ahlqvist**

[therese.ahlqvist@panthera.se](mailto:therese.ahlqvist@panthera.se)

Tel: +46 8-564 738 52