

S3 BRAKE FOR PUNCTURE FREE TYRE



S3 BRAKE



PUNCTURE FREE TYRE



BRAKE BLADE KIT

To improve the function of the S3 brake together with puncture free tyres (PU-tyre) we have developed a "kit", containing 2p brake blades made of steel and 6p binding screws, that is easily mounted on the S3 brake and increase the brake pressure considerable on the relatively "hard and slippery" puncture free tyre.

S3 Brake art.no. 4630000 Brake S3 cpl pair
PU Tyre art.no. 1000115 Tyre 24"x1" punct. free

Brake blade-kit art.no. 4630150 Brake plate S3 PU-Tyre

[Link to assembly video](#)