

UPGRADED S3 BACKREST

Click on images for larger format



Upgraded backrest with threaded knob



Reinforced backrest

UPGRADED S3 BACKREST


Effective from 2 October 2017, Panthera will be progressively upgrading the backrest on its S3 models..

This is to better secure the push handle attachment and reinforce the connection of the crossbar to the backrest tubing.

Please note that only the new version of the push handle and torso support attachment and can be mounted on this backrest. , see [info sheet 2017-14](#) and [info sheet 2017-18](#)

The price for the upgraded backrest will remain the same, but can be identified by the "ver.2" suffix on the label of the backrest frame and to the art.no.

Previous backrest



Previous backrest with unthreaded knobs



Previous crossbar