

U3 LIGHT BACKREST ON PANTHERA X



BACKREST U3 LIGHT



BACKREST PANTHERA X

USING THE U3 LIGHT BACKRESTS ON PANTHERA X

We would like to inform about the possibility to use the U3 Light backrests on Panthera X.

For those customers that wishes to use the backrests for U3 Light on Panthera X, it is possible to mount these, with a **smaller adaptation, made here at Panthera**. The only restriction is that the U3 Light backrest is not possible to mount in the rear-most backrest position, pos. 1, i.e. the position where the wheelchair is most easy balanced.

Valid for these backrests: (click on the images for more info)



U3 LIGHT STANDARD



FRAME-FOR FIXED BACKREST



BACKWARD CURVING



TAPERED



FLARED 3 CM



TAPERED 6 CM



POS. 1 NOT POSSIBLE TO USE

U3 LIGHT BACKREST ON PANTHERA X