

NEW MATERIAL IN THE OLD STANDARD REAR AXLE

Click on images for larger format



REAR AXLE IN ALUMINUM



REAR AXLE IN CARBON FIBRE

NEW MATERIAL IN THE OLD STANDARD REAR AXLE

We have updated the material of the tube on the old standard rear axle, with carbon fibre instead of aluminum that we used before.

This change is done continuously. First, the widths listed below are implemented, then the other widths will gradually be updated.

This change implicate that the older standard rear axle will have the same material as the present standard rear axle.

It gets lighter and stiffer.

The old standard rear axle is used together with the [assistant brake](#) and the [one arm drive](#), and some [power wheels](#).

Article numbers and price remains the same.

Valid for art.no:

UPDATED		
ART.NO.	DESCRIPTION	CHANGE
3301133	Tube rearaxle w bult 33	tube in carbon fiber instead of aluminum
3301136	Tube rearaxle w bult 36	tube in carbon fiber instead of aluminum
3301139	Tube rearaxle w bult 39	tube in carbon fiber instead of aluminum
3301145	Tube rearaxle w bult 45	tube in carbon fiber instead of aluminum