

New article from Panthera

Click on images for larger format



Knob anti-tip device S3- standard



New Knob anti-tip device S3- low profile



The knobs side by side

Knob anti-tip device S3 - low profile

We have designed a new knob for the [anti-tip device S3](#) that has a lower profile than the standard knob.

For wheelchairs with low chassis like [S3 short-low](#) and [S3 swing-low 0°](#) in combination with [backward curved back](#) rest the space around the anti-tip device knob is reduced. With this knob you can, in combinations above, still make use of the function of independently activate the anti-tip device sitting in the wheelchair.

When ordering anti-tip device S3 cpl with a low profile knob, we acknowledgement the order with article number for "Anti-tip device S3 pair" together with two "Knob anti-tip device S3 low profile". The knobs is delivered mounted on the anti-tip device and is included in the price for Anti-tip device complete.

The price for art.no. 4353121 is the same as for art.no. 4353122

[New article number](#)

art.no.	Description
4353121	Knob anti-tip device S3 low profile