


Bolt to inner casing castor fork X

Click on images for larger format



BOLT INNER CASING



BOLT INNER CASING

Bolt to inner casing castor fork on X

We have received comments from our customers that the bolt to the inner casing on the [Panthera X castor fork](#) can easily be damaged during service or repair. We therefore wish to make an improvement and replace this bolt with a new with a Torx T20 head.

The function and mounting is the same but with a reduced risk of damage the bolt when loosening and tightening.

We introduce this bolt when mounting on Panthera X wheelchairs and it is possible to order as a sparepart from 1st of October 2020.

The price remains the same.

Present bolt:

art.no. 6003408

Replacing bolt:

art.no. 6003409DL