

keeps you moving

## New wheelchair table - "Table S3"

Click on images for larger format



TABLE S3

### Now it's here!

**Finally we can offer a table for our S3 series! A flexible table for the active user.**

Can be used on all Panthera wheelchair models where S3 armrests can be mounted. One size table that fits seat width 33-50 cm.

Easy to mount on the armrests. The table can be mounted in 2 depths on the attachment sides. Height and depth position can also be further adjusted with the different positions of the armrests.

Stable construction that sits firmly in place ([Link to movie](#))

Maximum load 10 kg.

**Available for sale from today 13/12**

**Please contact Panthera if you have any questions**

Export / Marketing Manager

**Helena Lemoine**

Tel: +46-8-557 608 77

E-mail: [helena.lemoine@panthera.se](mailto:helena.lemoine@panthera.se)