

keeps you moving

## News! Footrest Carbon, for high mounting

Click on images for larger format



### Footrest Carbon, for high mounting

We have previously made adapted footrest in Carbon when the Bambino footrest has not been an alternative.

Now we can offer Footrest Carbon as an alternative in the standard range. It is made of a carbon plate, 140 mm deep, with two aluminum sides, mounted with clamps on the frame.

Available in three width, adapted designs for different chassis fronts; S3, U3 and X. For wheelchairs with abduction front you choose an S-type, two sizes larger than the seat width.

On X the Footrest Carbon needs to be use together with the standard footrest for stability.

Weight approx. 400 g

Description	Art.no.	Description	Art.no.
Footrest S3 carbon W. 27 cpl.	6504027	Footrest X carbon W. 33 cpl.	6506033
Footrest S3 carbon W. 30 cpl.	6504030	Footrest X carbon W. 36 cpl.	6506036
Footrest S3 carbon W. 33 cpl.	6504033	Footrest X carbon W. 39 cpl.	6506039
Footrest S3 carbon W. 36 cpl.	6504036	Footrest X carbon W. 42 cpl.	6506042
Footrest S3 carbon W. 39 cpl.	6504039	Footrest X carbon W. 45 cpl.	6506045
Footrest S3 carbon W. 42 cpl.	6504042		
Footrest U3 carbon W. 33 cpl.	6505033		
Footrest U3 carbon W. 36 cpl.	6505036		
Footrest U3 carbon W. 39 cpl.	6505039		
Footrest U3 carbon W. 42 cpl.	6505042		
Footrest U3 carbon W. 45 cpl.	6505045		

**Please contact Panthera if you have any questions**

Export / Marketing Manager

**Helena Lemoine**

Tel: +46-8-557 608 77

E-mail: [helena.lemoine@panthera.se](mailto:helena.lemoine@panthera.se)