

News - Cover plug Armrest S3

Click on images for larger format



REMOVE ARMREST S3



ARMREST S3 REMOVED



COVER PLUG ARMREST S3



COVER PLUG MOUNTED

Cover plug Armrest S3

Starting from the first of May 2022 Panthera now can offer the new product "Cover plug Armrest S3".

When the user needs to remove the armrest from the side guard there will be an opening in the side guard where the armrest was mounted.

With the Cover plug Armrest S3 you can now cover this opening and increase the comfort if there is a need to use the armrest as a sideguard only.

Art.no.	Description
4401140	Cover plug Armrest S3 pair
4401141	Cover plug Armrest S3 Left
4401142	Cover plug Armrest S3 Right

Please contact Panthera if you have any questions.

Export / Marketing Manager

Helena Lemoine

helena.lemoine@panthera.se

Tel: +46-8-557 608 77