

reTyre TRACTION 24" - now also available for 25" and 26" rear wheels



reTyre - TRACTION



reTyre - TRACTION

## World unique product, reTyre - TRACTION

Panthera promotes a world unique product, reTyre Traction tyre skins. **Now also available for 25" and 26"** rear wheels.

The reTyre Traction tyre skin is a detachable tyre cover that **gives better** friction and grip on soft / wet and slippery surfaces.

Designed for wheelchair driving with a pattern optimized for low rolling resistance , easy propulsion and braking.

The mounting is done easily with a velcro covered clip that is "clicked" on the rim, these are easy removable and mounted by the user.

The reTyre Traction tyre skins are designed to be used on manual wheelchairs with **24" , 25" and 26" rear wheels.**

The price is the same for 25" and 26" as for 24".

### [SEE MORE HERE](#)

Description	Art.no.
Tyre reTyre Traction 24" cpl pair	1100300
Tyre reTyre Traction 25" cpl pair	1100305
Tyre reTyre Traction 26" cpl pair	1100306

**Please contact Panthera if you have any questions.**

*Export / Marketing Manager*

**Helena Lemoine**

[helena.lemoine@permobil.com](mailto:helena.lemoine@permobil.com)

Tel: +46-8-557 608 77