



Document	Doc. Id	Version	Published Date	Page
CA & Compatibility Statement TMPL	CS-MK-PANTHERA	1	2023-07-07	1(4)
Issued by (Name)	Francesco Barbon	Date	2023-07-07	
Reviewed by (Name)	Mikael Beskow	Date	2023-07-07	
Approved by (Name)	Luciano Nosella	Date		

Compatibility Statement

MOTOTRONIK combined to PANTHERA wheelchairs



Permobil uses this compatibility statement (MDR art. §22) to list the wheelchairs to which the Mototronik can be mounted without compromising the integrity of the two devices and the safety of the user as long as the Mototronik is assembled by authorized and specialized personnel according to the instructions of the Mototronik Service Manual and as long as the user respects the instructions provided in the User's Manual of both devices.

Compatibility Statement between	
Party A	Rehateam s.r.l. Vicolo Negrelli 5, 31038 Castagnole di Paese TV - Italy
Party B	PANTHERA wheelchairs

Products of party A	
MOTOTRONIK 12" (including accessories)	MOTOTRONIK 14" (including accessories)
	
<p>The device is attached to the wheelchair by means of a connection frame. The clamps fixed to the wheelchair depend on the model of the same wheelchair.</p>	




Products of party B		for clamps mounting instructions consult the Mototronik Service Manual	
PANTHERA wheelchair (in alphabetical order)		PICTURE OF CLAMPS	PART NUMBER OF CLAMPS
S2			MO-AGL-DX-PSU MO-AGL-SX-PSU
S3			MO-AGL-DX-PSU MO-AGL-SX-PSU
U2			MO-AGL-DX-PSU MO-AGL-SX-PSU
U2 LIGHT			MO-AGL-DX-PSU MO-AGL-SX-PSU
U3			MO-AGL-DX-PSU MO-AGL-SX-PSU

U3 LIGHT			MO-AGL-DX-PSU MO-AGL-SX-PSU
----------	---	--	--------------------------------

Assembly instruction

For the models S2, S3, U2, U2 LIGHT, U3 and U3 LIGHT, the clamps are made of two parts to be fixed around the frame tube between the point **A** (bolt fixing the footplate tube) and the point **B** (beginning of the curve of the frame).



... linked document shows valida

Conclusion

The Mototronik and each product present in this compatibility list is in compliance with the MDR.

The use of product combinations is assessed to be within each manufacturer's respective product's intended purpose.

Record History

Compatibility Statement CS-MK-PANTHERA

Version	Change	Issued by
1	First issue	Francesco Barbon