

Adjusted by:

E-mail:

Adjusted for:

Telephone:

Date:

Save the guide to the desktop and open it from there.

Bambino 3 has an adjustable rear axle. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". Bambino 3 has a seat height of 41 cm in the front and 39 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the list of options.
Max userweight: 60 kg for seat width 24-33

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour

- Blue** - Racing Blue
- Black** - Anodic Black
- Clear

Optional colour - A higher price and

- longer delivery time for optional colour, Fill in art.no. in the square



2) Seat width (BB)

[art.no.](#)

- Clear

- 24 cm
- 27 cm
- 30 cm
- 33 cm

Art.no includes chassis, seat, backrest lock, anti-tip device



3) Backrest adjustable height

- 20-28 cm blue** 523BB17
- 20-28 cm black** 523BB17BL
- Clear
- 27-35 cm blue** 523BB18
- 27-35 cm black** 523BB18BL



4) Rear axle

- Carbon** 6 dgr 33030BB
- Rensa

- Carbon** 6 dgr +10 mm 33031BB

With 20" rearwheels +120 mm castors

- Other** rear axle or seat height

[Increase the seat height](#)

art. _____



5) Castor and fork

- 120 mm** [art.no.](#)
 - Clear
- With 20" rear wheels you need to raise the rear axle 10 mm

- 90 mm** [art.nr.](#)
- Standard with 20" rearwheels

- Other castor fork**

[Excel guide](#) art. _____



6) Footrest

- Footrest B3** - foldable, angle adjustable 42341BB
- Clear



7) Brake

- Brake B3** 4630200
- Clear

- Brake blade** B3 PU tyre SET 4630234



8) Rear wheel - without pushrim

- 22" Standard** [art.no.](#)
- Clear

- 20" Standard** [art.no.](#)
- 24" Standard** [art.no.](#)

[Rear wheel size](#)

- Other rear wheels**

art. _____



Tire - choose tire

- Schwalbe Right Run** 20",22",24"
- Clear

- Rugged** 24"
- Schwalbe **Marathon Plus** 22",24"

- PU** 20",22",24" (massive)



9) Pushrim

- Titanium** 190DD11 (x2 / wheelchair)

- Paragrip** 19051DDBT x2
- Tetragrip** 19052DDBTx2

- Clear DD = Wheel diameter

- Maxgrip 24"** (contains latex)
1903070 x2 (20")
1903080 x2 (22")
1903060 x2 (24")

Accessories

Side guard



- [Side guard B3](#), 4235000
- Clear

Arm support



- [Arm support B3](#) with side guard 4401400
- Adjustable in height in 3 pos.

Seat width adaptor



- [Seat width adaptor](#)-pair 4235020
 - Seat width adaptor-left 4235021
 - Seat width adaptor-right 4235022
 - Clear
- Limits the back rest angle forward to 16 mm. Narrows the seat width 7,5mm on each side.

Lateral support



- [Lateral support B3](#) adjustable cpl 4236000 /pcs
- Left Right
- Position
- Standard Standard
- Narrow -9mm Narrow -9mm

Push bar



- [Push bar B3](#) 4233000
- Height, angle adjustable with Quick release

Head support



- [Head support B3 Blue](#) 4710100
 - Head support B3 Black 4710100BL
 - Clear
 - Position Low Position High
- Height, angle and depth adjustable

Cushion Wedge cushion



- [2,5 cm](#) 572BB34
- [5 cm](#) 575BB34
- Clear
- [Wedge cushion](#) 57710BB not width 24 cm

Attachment pelvic belt

- Attachment plate pelvic belt cpl 9900247 backrest 2 point attachment
- Attachment clamp pelvic belt pair 9900248 2-point attachment
- [Attachment pelvic belt 4-p](#) plate + clamp 9900250 4-point attachment
- Attachment pelvic belt 4-p with clamps 9900248 x 2 4-point attachment

Spoke protectors



- Clear
- [24" Panthera logo](#) 4500000
- [22" Panthera logo](#) 4500100
- [20" Panthera logo](#) 4500021
- [24" Without logo](#) 4500011
- [22" Without logo](#) 4500111
- [24" Splash](#) 4500324
- [22" Splash](#) 4500322
- [20" Splash](#) 4500320
- [24" Car rim](#) 4500124
- [22" Car rim](#) 4500122
- [20" Car rim](#) 4500120
- [24" Animal Planet](#) 4500424
- [22" Animal Planet](#) 4500422
- [20" Animal Planet](#) 4500420
- [24" Smiley](#) 4500224
- [22" Smiley](#) 4500222
- [20" Smiley](#) 4500220

* Tube protection front; Can be too long if footrest is in its high position. Fits when rear axle is at std balance and footrest at low pos. 16 cm or lower, high pos. 28 cm or lower, from floor to top of footplate.

Other

- [Tube protection front](#) (chassis tube) * 9900251 (2 pcs)
 - [Tube protection rear](#) (backrest tube) 9900251 (1 pcs) Not with pushbar
 - [Backrest wedge](#) _____ pcs 4710012
 - [Tray-table](#) 6810000
 - [Magnetic holder](#) on rearwheel for tray table 6810010
 - [Torso strap](#) 4380238
- Chose art.no. for Freewheel according to the **footrest height** of the wheelchair
Read more about Freewheel on the [homepage](#)
- [Freewheel](#)
 - Freewheel-kit assembly parts 42342BB

Settings


Backrest angle



Note the measure (A) for desired backrest angle. 14-28mm. Standard=18 mm

A= _____ mm

Rear axle adjustment



Note the measure (A) for desired balance. Standard=77 mm

A= _____ mm

Footrest position



Note the height measure from the floor, (A) or (B) for the footrest.

[Bambino footplate with Feewheel - info](#)
A= _____ mm B= _____ mm

At 22" rear wheels
Measure high pos. 210-310 mm
Measure low pos. 150-220 mm

- Note depth-position: Std position Moved forward position 30 mm
- Note height-position: High position Low position

Transport marking

- YES do mark the wheelchair with attachment points for transport in disabled service Art.no.9000152

NOTE! See the manual for more information.

Extra