



















Adjusted by:

Tel:

E-mail:

Date:

| | | | |
|---|---|--|---|
| <p>Colour: <input type="radio"/> Clear <input type="radio"/> Red <input type="radio"/> Black</p> | <p>NOTE! Abduction front and seat width 24 cm only in red!</p> | | |
| <p>Seat width (BB): <input type="radio"/> Clear</p>  | <p><input type="radio"/> 24 cm * 356 24 00G <input type="radio"/> 24 cm Abd.* 357 24 00G</p> | <p><input type="radio"/> 27 cm 356 27 00(BL)G <input type="radio"/> 27 cm Abd.* 357 27 00G</p> | <p><input type="radio"/> 30 cm 356 30 00(BL)G <input type="radio"/> 30 cm Abd.* 357 30 00G</p> <p><input type="radio"/> 33 cm 356 33 00(BL)G <input type="radio"/> 33 cm Abd.* 357 33 00G</p> <p>NOTE! Use "BLG" after art.no. when choosing black colour</p> <p>* Only available in red</p> <p>Width 24 cm has the backrest moved 5 cm forward Width 24-33 cm max. user weight 60 kg art.no. includes: seat, backrest lock, foot rest, calf strap, brake and anti tip device.</p> |
| <p>Rear axle: <input type="radio"/> Clear</p>  | <p><input type="radio"/> Camber 6° 336 BB 16 BB is replaced by seat width, see "Seat width, above"</p> <p><input type="radio"/> other rear axle See list of options.</p> | <p>Rear axle adjustment balance of the chair</p> <p><input type="radio"/> std position, 22 mm <input type="radio"/> std position, 52 mm, backrest moved forward, width 24 <input type="radio"/> other position</p> <p>Note measure (A) for desired position. A= _____ mm</p> |  |
| <p>Castors: <input type="radio"/> Clear</p>  | <p><input type="radio"/> 120 mm for 20" rear wheels 210 00 25 <input type="radio"/> 120 mm for 24" rear wheels 210 01 05</p> | <p><input type="radio"/> 120 mm for 22" rear wheels 210 00 45</p> <p>*NOTE! You need two castors for each chair, the art.nr. is for one castor!</p> | |
| <p>Backrest height: <input type="radio"/> Clear</p>  | <p><input type="radio"/> Adjustable height 17-25 cm 323 BB 17</p> | <p><input type="radio"/> Adjustable height 27-35 cm 323 BB 18</p> | |
| <p>Footrest: <input type="radio"/> Clear</p>  | <p><input type="radio"/> Footplate, Bambino <input type="radio"/> Footplate swing up, Bambino 3 Read info sheet</p> <p>(40021BB) included in art.no. 35XBB00G, see "Seat width" (42341BB) Footplate Bambino 3 replace Footplate Bambino Swing Up in all widths except for Straight Front width 24 and 27cm you use: 4002325 width 24 cm and 4002328 width 27 cm</p> <p>When ordering Freewheel, always note the footrest height from the floor, in cm. <input type="text"/></p> | | |
| <p>Rear wheel: <input type="radio"/> Clear</p>  | <p>art.no.</p> <p><input type="radio"/> 20" Rear wheel Standard without pushrim</p> <p><input type="radio"/> 22" Rear wheel Standard without pushrim</p> <p><input type="radio"/> 24" Rear wheel Standard without pushrim</p> | <p>Pushrim: <input type="radio"/> Clear DD = wheel diameter</p> <p><input type="radio"/> Pushrim Titanium 190 DD 11</p> <p><input type="radio"/> Pushrim Paragrip 190 51 DDBT</p> <p><input type="radio"/> Pushrim Tetragrip 190 52 DDBT</p> <p><input type="radio"/> Pushrim Maxgrip* 190 30 (70) 20" (80) 22" (60) 24"</p> <p><input type="radio"/> other rear wheels <input type="text"/></p> | |
| <p>Tire: <input type="radio"/> Clear</p>  | <p><input type="radio"/> Tire - Right Run 20", 22", 24"</p> <p><input type="radio"/> Tire - Marathon Plus 22", 24"</p> <p><input type="radio"/> Tire - PU-massive 20", 22", 24"</p> | <p>NOTE! You need two rear wheels for each chair, the art.no. is for one rear wheel!</p> <p>*NOTE! Contains latex!!</p> | |
| <p>Spoke protectors: <input type="radio"/> Clear</p>  | <p><input type="radio"/> 24" Spoke protectors 450 00 00 with logo <input type="radio"/> 24" Spoke protectors 450 00 11 without logo</p> | <p><input type="radio"/> 22" Spoke protectors 450 01 00 with logo <input type="radio"/> 22" Spoke protectors 450 01 11 without logo</p> | <p><input type="radio"/> 20" Spoke protectors 450 00 21 with logo</p> |

| | |
|--|--|
| <p>Brakes:</p>  | <p><input type="radio"/> High brake, moved ahead</p> <p>(460 00 50) included in art.nr.356 BB 00G, see "Seat width"</p> |
| <p>Push handles:</p> <p><input type="radio"/> Clear</p>  | <p><input type="radio"/> Pushbar max height 100 cm with 22" wheels</p> <p>422 00 BB max. user weight 60 kg</p> <p><input type="radio"/> Height adjustable +10 cm</p> <p>425 10 00 max. user weigh 60 kg</p> <p><input type="radio"/> Height adjustable</p> <p>425 00 00 max. user weigh 60 kg</p> <p><input type="radio"/> Height adjustable +20 cm</p> <p>425 20 00 max. user weigh 60 kg</p> |
| <p>Anti tip device:</p>  | <p><input type="radio"/> Anti tip device, Bambino</p> <p>(435 10 01) included in art.no. 35X BB 00G, see "Seat width"</p> |
| <p>Sideguards:</p> <p><input type="radio"/> Clear</p>  | <p><input type="radio"/> Sideguards Bambino 20"</p> <p>441 00 15</p> <p><input type="radio"/> Sideguards Bambino 22"</p> <p>441 00 10</p> <p><input type="radio"/> Sideguards Bambino 24"</p> <p>441 04 10</p> <p><input type="radio"/> Sideguards Bambino winter 20"</p> <p>441 01 15</p> <p><input type="radio"/> Sideguards Bambino winter 22"</p> <p>441 01 10</p> <p><input type="radio"/> Sideguards Bambino winter 24"</p> <p>441 03 10</p> |
| <p>Cushions:</p> <p><input type="radio"/> Clear</p>  | <p><input type="radio"/> Cushion Micro long 2,5 cm for seat width 24 cm</p> <p>472 41 27</p> <p><input type="radio"/> Cushion Bambino 2,5 cm</p> <p>472 30 BB adjustable seatdepth 23-30 cm</p> <p><input type="radio"/> Cushion Micro long 5 cm for seat width 24 cm</p> <p>475 41 27</p> <p><input type="radio"/> Cushion Bambino 5 cm</p> <p>475 30 BB adjustable seatdepth 23-30 cm</p> |
| <p>Accessories:</p> | <p><input type="radio"/> Torso strap Bambino</p> <p>438 02 38</p> <p><input type="radio"/> Tube protection _____ pcs</p> <p>990 02 51 protects tube and paint</p> <p><input type="radio"/> Head cushion</p> <p>471 00 BB</p> <p><input type="radio"/> Wedge cushion</p> <p>577 10 BB</p>  |
| <p>Back rest angle:</p> <p>adjustment</p> | <p>Back rest angle adjusted with screw</p> <p>Note measure (A) for desired angle A= _____ mm</p>  |
| <p>Footrest:</p> <p>adjustment</p> | <p>Footrest adjustments</p> <p>Specify measures (A and B) for the footrest A= _____ B= _____ mm</p> <p>Standard measures A=141 mm and B=120 mm at 20" rear wheels Standard measures A=165 mm and B=150 mm at 22" rear wheels Standard measures A=190 mm and B=175 mm at 24" rear wheels</p>  |
| <p>Marking:</p> | <p><input type="radio"/> YES mark the wheelchair for transportation in disabled service vehicles.</p> <p><input type="radio"/> NO do NOT mark the wheelchair for transportation in disabled service vehicles</p> <div style="border: 1px solid black; padding: 2px; display: inline-block;"> <p>NB! See instruction for use.</p> </div> <p>Art.no.9000152</p>  |
| <p>Other:</p> | |