

Utprovad av:

E-post:

Utprovad för:

Telefon:

Datum:

Total weight (g):

The weight changes based on the choices made below.

S3 short Abd has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 Short Abd has a seat height of 45 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.. **S3 short Abd is 50 mm shorter than S3 and has an abducted front and 4° seat angle.**
Max userweight is 100 kg for seat width 27-33 cm.

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function
Design your wheelchair by clicking in the circles in this guide.
Blue, underlined text is clickable to **show info** or image.
BB = seatwidth, **HH**= backrest height

Original design

Alternative options

Other options

[info](#)

1) **Chassis-colour** [Click on the link to see up-to-date colours](#)



2) **Seat width (BB)** 27 cm 30 cm 33 cm
[art.no.](#)
 Clear

Art.no include: chassis, seat, calf strap,

3) **Seat depth** Adjustable seat depth 27-33 cm
 Clear



4) **Back rest** [Height \(HH\)](#)
 [Standard](#) 523BBHH VER2 35 cm
 25 cm 40 cm
 30 cm
 Clear

[Alterantive, back rests art.nr.](#)
art. _____



5) **Rear axle**
 [Bambino 6 dgr](#) 336BB16
 Clear

[Carbon](#) +15mm 33023BC* * seat width **BB** +1

[Other](#) rear axle or seat height
[Increase the seat height](#)



6) **Castor and fork**
 [120 mm art.no.](#)
 Clear

[90 mm art.no.](#)
 [150 mm art.no.](#)



7) **Footrest**
 [Footrest Titanium](#), small heel bar
 Clear

[Bambino footplate art.no.](#)
 [Foldable footplates art.no.](#) Wheelchair with abduction front (Abd): Choose footplate 2 sizes larger than the seat width.
 [Heel strap for fotplates art.no.](#)

[Footrest](#) extended reinforced [art.no.](#) [Needed for Freewheel](#)



8) **Brake**
 [Brake S3](#) 4630000
 Clear

[High brake, moved ahead](#) 4600050

[Assistant brake](#)
 [Brake S3](#) moved rearward 4630124
 [Brake S3](#) moved forward 4630010



9) **Rear wheel - without pushrim**
 24" [standard art.no.](#)
 Clear

22" [standard art.no.](#)
Rear wheel sizes

[Other rear wheel](#)



Tire- choose tire
 [Schwalbe Right Run](#) 22", 24" Standard
 Clear

[Rugged](#) 24"
 Schwalbe [Marathon Plus](#) 22",24"

[PU](#) 24" (massive)
 [PU](#) 22" (massive)
 [reTyre](#)-tyre skins (does not fit PU tyres)



10) **Push rim**
 [Titanium](#) 190DD11 (x2 / wheelchair)
 Clear **DD** = Wheel diameter

[Aluminum X, 24"](#) 1904404 x2
 [Paragrip](#) 19051DDBT x2
 [Tetragrip](#) 19052DDBT x2 **DD** = Wheel diameter

[Maxgrip](#) 24" (contains latex) 1903060 x2
 [Maxgrip](#) 22" (contains latex) 1903080 x2

Accessories



Side guard

- [Side guard S3](#) [Video](#)
art.no.
- Clear



Arm support

- [Arm support S3](#) [Video](#)
art.no.



Side guard Carbon

- [Carbon](#)
art.no.
- [Carbon , small cover](#)
art.no.



Push handles H/S

- [Height adjustable S3](#) [Video](#)
4252060 VER2
 - [Height adjustable S3+10 cm](#)
4252063 VER2
 - Clear
- Max. userweight 100 kg



Push bar

- [Push bar S3](#)
42400BB VER2



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Anti tip

- [Anti tip](#) std seat height
4353000 [Video](#)

[Other seat height](#)



Cushion

- | | |
|---|---|
| Short | Short adjustable BB 30-33 |
| <input type="radio"/> 2,5 cm
572BB35 | <input type="radio"/> 2,5 cm
47234BB |
| <input type="radio"/> 5 cm
575BB35 | <input type="radio"/> 5 cm
47534BB |
| <input type="radio"/> Clear | |



Spoke protectors

- | | |
|--|--|
| <input type="radio"/> 24" logo
4500000 | <input type="radio"/> 22" logo
4500100 |
| <input type="radio"/> 24" no logo
4500011 | <input type="radio"/> 22" no logo
4500111 |
| <input type="radio"/> Clear | |




Other

- [Tube protection](#) _____ P
9900251 (x3/wheelchair)
- [Crutch holder S3](#) _____ P
9900301 (x1/krycka) [Video](#)
- [Backrest wedge](#) _____ P
4710012

Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)

- [Freewheel](#)



Wedge cushion

- [Wedge cushion](#)
57710BB

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Std with 6° rear axle =31 mm,
Std with Carbon rear axle = 40 mm



A= _____ mm



Footrest height [Footbow](#) - [Footplates](#) - [Bambino plate](#)

Note measure (A) and (B) for desired footrest .

A= _____ mm B= _____ mm



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle

Art.no.9000152

NOTE! See [manual](#) for more information.

Extra