

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

Vikten ändrar sig beroende på vilka val du gör nedan.

S3 Swing +25 has divided, swing away leg support and an adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 Swing extended 25 is 25 mm longer than S3 Swing and has a seat height of 45cm in front and 43cm at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.

Max userweight: 100 kg for seat width 36-42 and 150kg for seat width 45-50 cm.

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. **BB** = seatwidth, **HH**= backrest height

Original design

Alternative options

Other options

[info](#)

1) [Chassis-colour](#) [Click on the link to see up-to-date colours](#)



2) [Seat width \(BB\)](#)

[art.no.](#)

- 36 cm 39 cm 42 cm 45 cm 50 cm

Art.no include: chassis, seat, backrest lock

Clear

3) [Seat depth](#)

[art.no.](#)

- 35 cm 40 cm

Clear

- 37,5 cm 42,5 cm

4) [Back rest](#)

[Standard](#)

523BBHH VER2

Clear

Height (HH)

20 cm 35 cm

25 cm 40 cm

30 cm 45 cm

Clear

[Tapered](#)

525BBHH VER2

[Curved back](#)

524BBHH VER2

[Info](#)

[Back rest, extended range](#)

[Frame fo fixed backrest systems](#)

art. _____

5) [Rear axle](#)

[Carbon standard](#)

33023BB

Clear

[Carbon +15 mm](#)

33023BC*

* seat width **BB** +1

Use this rear axle when chosing arm support at Accessories and/or increased seat height

[Other](#) rear axle or seat height

[Increase seat height](#)

6) [Castor and fork](#)

[120 mm](#)

Clear

[90 mm](#)

[150 mm](#)

7) [Footrest](#)

[Swing S3 cpl pair](#)

[art.no.](#)

Clear

[Footplate large](#)

[art.no.](#)

[Extended tubes 40 mm](#)

[art.no.](#)

[Heel straps for footplates](#)

[art.no.](#)

8) [Brake](#)

[Brake S3](#)

4630000

Clear

[Other brake](#)

9) [Rear wheel - without pushrim](#)

24" [Standard](#)

[art.no.](#)

Clear

22" [Standard](#)

[Rear wheel sizes art.no.](#)

25" [Spinergy X](#) black spokes

26" [Spinergy X](#) black spokes

[Other rear wheel](#)

Tire - choose tire

[Right Run](#) 22",24",25",26"

Standard

Clear

[Rugged](#) 24"

[Marathon Plus](#) 22",24",25",26"

[PU](#) 24" (massive)

[PU](#) 22" (massive)

[reTyre](#)-tyre skins (does not fit PU tyres)

10) [Pushrim](#)

[Titanium](#)

190DD11 (x2)

[Rensa](#)

DD = Hjulldiameter

[Paragrip](#)

19051DDBT x2

[Tetragrip](#)

19052DDBTx2

DD = Hjulldiameter

[Maxgrip](#) 24" (contains latex)

1903060 x2

[Maxgrepp](#) 22" (contains latex)

1903080 x2

Accessories



Side guard S3

- [Side guard S3](#) [Video](#)
[art.no.](#)
- Clear



Arm support

- [Arm support S3](#) [Video](#)
[art.no.](#)



Side guard Carbon

- [Carbon](#)
[art.no.](#)
- [Carbon, small cover](#)
[art.no.](#)
- [Carbon, large cover *](#)
[art.no.](#)
- * Only with one arm brake



Push handles H/S

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3+10 cm](#)
4252063 VER2
Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push bar

- [Push bar S3](#)
42400BB VER2

Max. userweight 100kg.



Tippskydd

- [Anti-tip](#) std seat height
4353000 [Video](#)

[Other seat height](#)



Cushion

- Seat depth 35 cm Seat depth 40 cm
- [2,5 cm](#) [2,5 cm](#)
572BB35 572BB40
- [5 cm](#) [5 cm](#)
575BB35 575BB40
- Clear




Spoke protectors

- [24" logo](#) [22" logo](#)
4500000 4500100
- [24" no logo](#) [22" no logo](#)
4500011 4500111
- Clear



Other

- [Tube protection](#) _____ p
- 9900251 (x3/wheelchair)
- [Crutch holder S3](#) _____ p
- 9900301 (x1/crutch) [Video](#)
- [Back rest wedge](#) _____ p
- 4710012



Wedge cushion

- [Wedge cushion](#)
57710BB


Set-up adjustments



Backrest angle

Ange måttet (A) för önskad ryggvinkel.
14-28 mm. Standard=18 mm

A= _____ mm



Rear axle adjustment

Ange måttet (A) för önskad balansering.
Standard=40 mm

A= _____ mm



Footrest height

Note measure (A) and (B) for the footrest.

A= _____ mm B= _____ mm

Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle.

Art.no.9000152

NOTE! See [manual](#) for more information.



Extra