

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (lb):

The weight changes based on the choices made below.

Panthera X has the chassis and backrest made in carbon fibre and a rigid rear axle. The balance is adjusted by positioning the backrest in 4 different positions. Panthera X is delivered without anti-tip device and is intended for advanced users only. Max user weight is 220 lb for seat width 13-18" cm.

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)



1) Seat width (BB) Chassis [art.no.](#)

- 13"
- 14"
- 15"
- 16,5"
- 18"
- Clear



2) Backrest height 8,7-13,8" (HH)

- 8,7"
- 9,05"
- 9,5"
- 9,8"
- 10,2"
- 10,6"
- 11"
- 11,4"
- 11,8"
- 12,2"
- 12,6"
- 13"
- 13,04"
- 13,08"
- Clear



3) Balance [info](#)

- Position 1**
Most easy balanced
- Position 2**
Less easy balanced
- Position 3**
Medium balanced
- Position 4**
Least easy balanced
- Clear

The positioning of the backrest on the chassis affects the balance / COG of the wheelchair. "easy balanced" = easy to practice rear wheel balance. This can later easily be changed.



4) Castor

- 3" for 24" rear wheel
2305130 (x2 / wheelchair)
- Clear

- 3" for 25" rear wheels
2305131 x2
- 3" for 26" rear wheels
2305132 x2

- 5" for 24" rear wheels
2305133 x2



5) Footrest

- [Footrest Alu](#)
36004 BB
- Clear

- [Footrest Titanium](#)
36915BB
- [Footrest Titanium, extended, reinforced](#)
36006BB [needed for Freewheel](#)
- [Plastic plate for foot bow](#)
36010BB



7) Brake [art.no.](#)

- [One hand brake, right](#)
- Clear

- [One hand brake, left](#)

- [Other brake](#)

art. _____



8) Rear wheel - without pushrim

- [24" Spinergy X](#) black spokes
- [24" Spinergy X](#) white spokes
- Clear

[art.no.](#)

- [25" Spinergy X](#) black spokes
- [25" Spinergy X](#) white spokes
- [26" Spinergy X](#) black spokes
- [26" Spinergy X](#) white spokes

[Rear wheel sizes](#)

- [Other rear wheel](#)

[List of options](#) art. _____



9) Tire

- [Schwalbe One](#) 24"
- Clear

- [Right Run](#) 24",25",26"
- [Marathon Plus](#) 24",25",26"

- [Rugged](#) 24"

[info](#)



10) Push rim

- [Aluminum X](#)
1904404 (5) (x2 / wheelchair) 24" (25")
- Clear

- [Titanium](#) **DD** = Wheel diameter in inches
190DD11 x2
- [Paragrip](#)
19051DDBTx2
- [Tetragrip](#)
19052DDBTx2

- [Maxgrip](#) 24" (contains latex)
1903060 x2

Accessories



Side guard

- [U2 light](#)
art.no. [info](#)
- Clear



Side guard

- [Carbon, small cover](#)
art.no.



Side guard

- [Carbon, large cover](#)
art.no.



Cushion

- X 2"** seat depth 15,8"
47501BB
- Rensa



Arm hook

- [Arm hook X](#)
4210100
- Rensa



Transport protection

- [Housing](#)
9024009
- [Rear end](#)
9024012
- [Rear axle](#)
9024008



Other

- [Freewheel](#) [info](#)
2600102
- [Tube protection](#)
9900252 (x3 / wheelchair) ____ p


Set-up adjustments



Backrest angle

Note measure (A) for desired back rest angle
0,39-0,78". Standard=0,7"

A= _____ "



Footrest height

Note measure (A) for desired footrest height
Standard = 2,75"

A= _____ "

Extra